

DENAS-PCM Pro

ILLUSTRATED ATLAS OF APPLICATION ZONES

English Companion to Part II - Instructions for Use

Source. The anatomical illustrations are extracted from Appendix B of the manufacturer's Russian DENAS-PCM Pro Operating Manual, TRTK 23.0-03.70-02 RE. Russian page text has been omitted and replaced here with English headings and captions.

Important. This is an English translation/companion document, not an official manufacturer-issued English IFU. Use only on intact skin and follow the safety instructions and contraindications in the Operating Manual.

Allergy

Zone 1: epigastric abdominal area. Zone 2: area of the 7th cervical vertebra.




Pain

Treat the painful area. For a large painful area, reposition the electrode sequentially.



Severe Pain

Treat the area of maximum pain.



VSD

Cervical-collar zone.



Bloating

Anterior abdominal wall.




Gynecologic Pain

Painful area of the lower abdomen.



Gynecology

Zone 1: soles. Zone 2: cervical-collar zone. Zone 3: lumbosacral zone.



Hypertension

Zone 1: area of C2. Zone 2: area of C7.



Hypotension

Cervical-collar zone.



Head

Zone 1: C2 area. Zone 2: cervical-collar zone. Zone 3: hand projection of the head according to Su-Jok.



Throat

Zone 1: hand projection of the neck according to Su-Jok. Zone 2: submandibular area.



Diarrhea

Zone 1: anterior abdominal wall. Zone 2: lumbosacral zone.



GI Pain

Painful abdominal area.



Constipation

Zone 1: anterior abdominal wall. Zone 2: lumbosacral zone.




Cough

Zone 1: submandibular area. Zone 2: jugular notch.



Muscles

Area of muscle pain or injury.





Muscles / Sports

Affected muscle area. See the next page for self-adhesive electrode placement examples.



Muscles / Sports - Electrode Placement

Examples of self-adhesive electrode placement on upper arm, trunk/hip, thigh, and lower leg.



Runny Nose

Zone 1: hand projection of the nose according to Su-Jok. Zone 2: nasal wings and projections of the paranasal sinuses.



Potency

Zone 1: lumbosacral area. Zone 2: cervical-collar area. Zone 3: feet.




Kidneys

Zone 1: lumbosacral area. Zone 2: feet.



Back

Cervical spine for neck pain; thoracic/lumbar/sacral spine for back pain.



Joints

Zone 1: affected joint. Zone 2: lower back for lower-limb joints or cervical-collar area for upper-limb joints.



Nausea

Epigastric area and right/left hypochondria.



Trauma

Treat around the affected area. Do not place electrodes directly over wounds or broken skin.



Choking / Dyspnea

Zone 1: jugular notch. Zone 2: area of C7.



Cellulite

Affected area; the source manual describes course treatment by individual body zone.



Cellulite - Electrode Placement

Examples of self-adhesive electrode placement for the abdomen, buttocks, outer thighs, posterior thighs, and anterior thighs.

